

## Neonatal Noninvasive Bilirubin Meter

## The BiliCare™ Transcutaneous Bilirubin Meter

The BiliCare System measures the level of bilirubin using a patented technology to provide accurate results. It's technology is based on spectroscopy.

The System measures the level of transcutaneous bilirubin in the newborn by transmitting light at different wavelengths through the outer ear. The amount of light absorbed by the bilirubin is calculated according to a customized algorithm.

- Less sensitive to motion artifacts
- Not sensitive to differences in pressure applied or angle of device

# The LED advantage

LEDs do not deteriorate over time. No shift in power or in wavelength means no need for routine calibration.

Optional calibration check allows user to confirm calibration when desired.

The calibration check process is quick and easy & provides immediate results.

## **Communication features**


BiliCare includes communication features that may connect it to a PC and/or to a HIS via HL-7

## **BiliCare Viewer**

Unique S/W that displays the results in tables and graphs.

These may be printed and exported.

Search by Patient ID1:		ō [	3 🖳	As .	0	
Potent RI - Sales	III. Section 200	Sicin Dates - Deliver 15 at last		SIALAR	San Marin	SO SO I
(8256 EZ)	Patient (DJ)	21: Own 13 18:89	1.5 mg/d.	o gigh respective -	SECRY SHE	348279
Miron III.	Avoies (02:	30-Dec-14-17-00	K.Tmgrds	Khile-Stills Percentifie	Micey S.W.	500269
And head \$250. Toronto.		Martin Daller of Switch and Sil	61.30	W (A) (A) (B)		
Marini No.	(National 2011)	29-Dati-19 (8-40)	HARWS.	Philly This Personality	BRCay NA	340219
Marke St.	Patient (Gd)	30 Dec 15 JE 18	7.2 mg/st.	Section 75th Percentile.	BBC219 1/70	341275
Name In	Patient Kids	49-Geo-15-14-15	254906	20th-00th Percentiles	Miller 5-76	148219
Political River Communication	1000000	Birth Defet 17 Sec-15 (0.00)	20 40	IN ALABA	TO THE REAL PROPERTY.	7707
Metar XX.	Pyleet \$10:	28-0ec-14-80.00	6.8 mg/st.	#-3586 Procuestie	MiCare Silv.	383399
there are	Patient Sig:	20 Dec 13-93-00	17.0 mg/dL	FOOT TOTAL PROCESSES	BRCare 3/RI	349238
Nova 30:	Patient IO2:	39-Dec-19 13:49	13.0 mg/d.	1900 Sonoelle	BBCon S.W.	148299
Marine ID:	Patient IDS:	48-Dec/ 15 19-00	Edward.	75th 15th Percentile.	MiCey SR	344219
Political III In Holyto	1	MICOLDINAL LIPERS-25 (ACM)	前 3	100 (A. (A. (A.))		
Horse Etc.	Patient EU:	13 Dec-10-18-03	55 mg/st.	F. TOTA FACINGIA	BROWN 5-76	341275
Marcon 80:	Patient Eld:	99-Dec-16-16-68	6.9 mg/d.	SRIDY 75th Percentille	BBCorn S/Ri	598279
Notice ID-	Pytieri ID2	12-0+c-12-22-24	0.2 mg/d.	th 15th Parcentlle	BBCare S.R.	045276





## www.mennenmedical.com

United States: Mennen Medical LLC. 868 W. Street Rd.

Unit #405, Warminster, PA 18974, USA Phone: +215 259-1020

Fax: +215 357-2010

Distributed in the UK by KONTRON

United Kingdom: Charter Kontron Ltd. Scorpio House, Rockingham Drive, Linford Wood, Milton Keynes MK14 6LY UK Phone: +44 (0)1908 646070 Fax +44 (0)1908 646030







**Neonatal Noninvasive Bilirubin Meter** 



Accurately and conveniently measures and monitors the level of bilirubin to assess the risk of neonatal hyperbilirubinemia.







Israel: Mennen Medical Ltd.

4 Paran Street, Building #11

Yavne, 8122503, Israel

Phone: +972-8-9323333

Fax: +972-8-9328510

## THE NEXT GENERATION TRANSCUTANEOUS BILIRUBIN METER

#### Simple for the User & Gentle on the Baby

- Reduces the need for Total Blood Bilirubin and the number of sticks in the baby's heel, by providing an accurate bedside measurement
- Only requires a single measurement per test

## LED Transmission Technology Reduces User Variances & Maintenance Requirements

- Less sensitive to motion artifacts & measurement differences based on user technique
- LEDs do not deteriorate over time eliminating the requirement for routine device calibration

## ADVANCED TECHNOLOGY FOR ACCURATE MEASUREMENTS

## Speed and convenience

- Reduces number of blood tests via heel puncture
- Can use just one measurement or average 2–3 measurements
- Can be used with or without entering baby's ID
- Large & clear display for easy use
- Resistive touch screen can be used with or without gloves
- · Menu-driven user interface minimizes training needs, facilitating use by multiple caregivers
- Memory up to 100 measurements saves time in transcribing and comparing results
- · Option to enter both patient and user identification facilitates hospital audits and quality assurance
- Barcode scanner for quick and accurate entry of caregiver and baby identification

#### Accuracy

- The BiliCare System was validated in clinical studies to produce results comparable to blood tests
- Repeatability tests confirm reliable results

## Low cost of ownership

- The BiliCare System is capable of working with or without disposable tip cover
- LEDs do not require routine calibration, minimizing maintenance and cost associated with service
- · Long-lasting, rechargeable battery with cradle allows extended use and ease of recharging



**Tip Covers** 


## **MEASURING WITH BILICARE**

Neonatal jaundice is a common condition, but when not carefully monitored, hyperbilirubinemia can lead to kernicterus, a devastating condition.

Early measurement of bilirubin levels, leading to early intervention is recommended but the time to take a blood sample and the trauma to the baby made universal screening difficult. Now, the BiliCare System provides a cost and time effective way to monitor bilirubin levels without unnecessary heel sticks.

Ease of Use - just follow the simple prompts for quick measurement where and when you need it:

Simply press measure on the home screen to initiate the measurement process



Position the BiliCare clip on the baby's ear to take your measurement



You are now ready to review your results or take additional measurements per your hospital protocol and device set-up

