



## Global Neonatal Care Expertise

• Beijing • Shenzhen • Hong Kong • Brisbane



## Conlett Pty Ltd

Unit 25, 17 Cairns St, Loganholme QLD Australia, 4129  
Tel: + 617 32099957 Fax : +617 3209976  
Mobile: + 61 413308385  
www.conlett.com

**Neo-Guard** Family is designated to protect sensitive neonatal nasal and facial skin during respiratory therapies, effectively guarding neonatal septum from trauma and irritation caused by pressure & motion of prongs and masks.



- Hydrocolloid material is super skin friendly even to the sensitive skin of neonates.
- Repositioning of the Neo-Guard on neonates is easy and does not cause any disconformity.
- The special hook & loop design at the bottom of every piece of Neo-Guard provides an extra method of securing multiple types and sizes of nasal device.
- Easy observation of the condition of the nasal nares.
- Multiple sizes & designs available for the best matchup of therapies and baby weighs.

**Optimization** Three unique designs best tuned for their designated oxygenation scenario. Pick the best match for the best result.



### Neo-Guard Plus

- Neo-Guard Plus integrates an extra layer of foam on top of the regular Neo-Guard design to further protect the nasal septum

### Neo-Guard

- Regular Neo-Guard is specially designed for nasal CPAP.
- It adopts a partially cut-through Star around the holes to generate retention on the prongs, effectively securing prong position and providing better seal for a stable alveolar pressure.
- Neo-Guard effectively holds the prongs by its retention and thus prevents intact between prongs and baby's nose



### High Flow Neo-Guard

- High Flow Neo-Guard is set to provide a clear flow path for high flow nasal cannula therapy, by having hole size much larger than that of the prongs.

## Ordering Information

NG000	CPAP Neo-Guard<700g
NG001	CPAP Neo-Guard 700g~1250g
NG002	CPAP Neo-Guard 1250g~2000g
NG003	CPAP Neo-Guard 2000g~3000g
NG004	CPAP Neo-Guard >3000g

NGP-S-M	CPAP Neo-Guard w/ foam<2000g
NGP-M-L	CPAP Neo-Guard w/ foam>2000g

HFNG-S	High Flow Neo-Guard <1250g
HFNG-M	High Flow Neo-Guard 1250g~2000Gg
HFNG-L	High Flow Neo-Guard>2000g